

# Human ATM/Serine-protein kinase ATM ELISA Kit

Catalog No: #EK19603



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ATM/Serine-protein kinase ATM ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Serine-protein kinase ATM
Target Name	Serine-protein kinase ATM
Other Names	Ataxia telangiectasia mutated,A-T mutated
Accession No.	Q13315

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.16ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ATM

Note: This product is for in vitro research use only and is not intended for use in humans or animals.