

Human BRCA2/Breast cancer type 2 susceptibility protein ELISA Kit

Catalog No: #EK19605



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK19605-1 48T #EK19605-2 96T

Description

Product Name	Human BRCA2/Breast cancer type 2 susceptibility protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Breast cancer type 2 susceptibility protein
Target Name	Breast cancer type 2 susceptibility protein
Other Names	Fanconi anemia group D1 protein,FACD,FANCD1
Accession No.	P51587

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40.6pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BRCA2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.