

Mouse Cdh5/Cadherin-5 ELISA Kit

Catalog No: #EK19606



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Cdh5/Cadherin-5 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Cadherin-5 |
| Target Name | Cadherin-5 |
| Other Names | Vascular endothelial cadherin, VE-cadherin, CD144 |
| Accession No. | P55284 |

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.33ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdh5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.