

# Rat Ada/Adenosine deaminase ELISA Kit

Catalog No: #EK19608



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Ada/Adenosine deaminase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Adenosine deaminase
Target Name	Adenosine deaminase
Other Names	Adenosine aminohydrolase
Accession No.	Q920P6

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.176ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Ada

Note: This product is for in vitro research use only and is not intended for use in humans or animals.