

Human CSN3/Kappa-casein ELISA Kit

Catalog No: #EK19613



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CSN3/Kappa-casein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Kappa-casein
Target Name	Kappa-casein
Other Names	CASK,CSN10,CSNK
Accession No.	P07498

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.066ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CSN3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.