

Human COQ2/4-hydroxybenzoate polyprenyltransferase, mitochondrial ELISA Kit

Catalog No: #EK19616

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human COQ2/4-hydroxybenzoate polyprenyltransferase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human 4-hydroxybenzoate polyprenyltransferase, mitochondrial
Target Name	4-hydroxybenzoate polyprenyltransferase, mitochondrial
Other Names	4-HB polyprenyltransferase,4-hydroxybenzoate decaprenyltransferase,COQ2 homolog,hCOQ2,Para-hydroxybenzoate--polyprenyltransferase,PHB:PPT,CL640
Accession No.	Q96H96

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.188ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:COQ2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.