

Mouse Chit1/Chitotriosidase-1 ELISA Kit

Catalog No: #EK19617



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Chit1/Chitotriosidase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Chitotriosidase-1
Target Name	Chitotriosidase-1
Other Names	Chitinase-1
Accession No.	Q9D7Q1

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.099ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Chit1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.