

Mouse Col4a3bp/Collagen type IV alpha-3-binding protein ELISA Kit

Catalog No: #EK19618



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Col4a3bp/Collagen type IV alpha-3-binding protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen type IV alpha-3-binding protein
Target Name	Collagen type IV alpha-3-binding protein
Other Names	Ceramide transfer protein, Goodpasture antigen-binding protein, GPBP, START domain-containing protein 11, StARD11, StAR-related lipid transfer protein 11, Stard11
Accession No.	Q9EQG9

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Col4a3bp

Note: This product is for in vitro research use only and is not intended for use in humans or animals.