

# Mouse Bcap31/B-cell receptor-associated protein 31 ELISA Kit

Catalog No: #EK19622



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bcap31/B-cell receptor-associated protein 31 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse B-cell receptor-associated protein 31
Target Name	B-cell receptor-associated protein 31
Other Names	BCR-associated protein 31,p28,Bap31
Accession No.	Q61335

## Application Details

Assay Type: Sandwich  
Detect Range:78-5000pg/mL  
Sensitivity:39.8pg/mL  
Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul  
Assay Time:3-5h  
Detection wavelength:450 nm  
gene\_full\_name:Bcap31

Note: This product is for in vitro research use only and is not intended for use in humans or animals.