

Mouse Clec3b/Tetranectin ELISA Kit

Catalog No: #EK19623



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Clec3b/Tetranectin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Tetranectin
Target Name	Tetranectin
Other Names	TN,C-type lectin domain family 3 member B,Plasminogen kringle 4-binding protein,Tna
Accession No.	P43025

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.152ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Clec3b

Note: This product is for in vitro research use only and is not intended for use in humans or animals.