

Human CST9/Cystatin-9 ELISA Kit

Catalog No: #EK19624



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CST9/Cystatin-9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cystatin-9
Target Name	Cystatin-9
Other Names	Cystatin-like molecule,CLM,CTES7A
Accession No.	Q5W186

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CST9

Note: This product is for in vitro research use only and is not intended for use in humans or animals.