

Human AKT1S1/Proline-rich AKT1 substrate 1 ELISA Kit

Catalog No: #EK19630



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human AKT1S1/Proline-rich AKT1 substrate 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Proline-rich AKT1 substrate 1
Target Name	Proline-rich AKT1 substrate 1
Other Names	40 kDa proline-rich AKT substrate, PRAS40
Accession No.	Q96B36

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AKT1S1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.