

Rat Aktip/AKT-interacting protein ELISA Kit

Catalog No: #EK19632



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Aktip/AKT-interacting protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat AKT-interacting protein
Target Name	AKT-interacting protein
Other Names	Fused toes protein homolog,Fts
Accession No.	Q5FVH4

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.085ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aktip

Note: This product is for in vitro research use only and is not intended for use in humans or animals.