

# Human AICDA/Single-stranded DNA cytosine deaminase ELISA Kit

Catalog No: #EK19633



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human AICDA/Single-stranded DNA cytosine deaminase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Single-stranded DNA cytosine deaminase
Target Name	Single-stranded DNA cytosine deaminase
Other Names	Activation-induced cytidine deaminase,AID,Cytidine aminohydrolase,AID
Accession No.	Q9GZX7

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.05ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:AICDA

Note: This product is for in vitro research use only and is not intended for use in humans or animals.