

# Mouse Agmat/Agmatinase, mitochondrial ELISA Kit

Catalog No: #EK19636



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Agmat/Agmatinase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Agmatinase, mitochondrial
Target Name	Agmatinase, mitochondrial
Other Names	Agmatine ureohydrolase,AUH
Accession No.	A2AS89

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Agmat

Note: This product is for in vitro research use only and is not intended for use in humans or animals.