

Human ACO1/Cytoplasmic aconitate hydratase ELISA Kit

Catalog No: #EK19637

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ACO1/Cytoplasmic aconitate hydratase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cytoplasmic aconitate hydratase
Target Name	Cytoplasmic aconitate hydratase
Other Names	Aconitase,Citrate hydro-lyase,Ferritin repressor protein,Iron regulatory protein 1,IRP1,Iron-responsive element-binding protein 1,IRE-BP 1,IREB1
Accession No.	P21399

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ACO1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.