

Rat Aco1/Cytoplasmic aconitate hydratase ELISA Kit

Catalog No: #EK19638

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Aco1/Cytoplasmic aconitate hydratase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cytoplasmic aconitate hydratase
Target Name	Cytoplasmic aconitate hydratase
Other Names	Aconitase,Citrate hydro-lyase,Iron regulatory protein 1,IRP1,Iron-responsive element-binding protein 1,IRE-BP 1,Ireb1,Irebp
Accession No.	Q63270

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.178ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aco1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.