

# Mouse Ca9/Carbonic anhydrase 9 ELISA Kit

Catalog No: #EK19643



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ca9/Carbonic anhydrase 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carbonic anhydrase 9
Target Name	Carbonic anhydrase 9
Other Names	Carbonate dehydratase IX,Carbonic anhydrase IX,CA-IX,Membrane antigen MN homolog,Car9
Accession No.	Q8VHB5

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.14ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Ca9

Note: This product is for in vitro research use only and is not intended for use in humans or animals.