

# Rat Arhgef7/Rho guanine nucleotide exchange factor 7 ELISA Kit

Catalog No: #EK19645



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Rat Arhgef7/Rho guanine nucleotide exchange factor 7 ELISA Kit       |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Rt   |
| Specificity        | Natural and recombinant Rat Rho guanine nucleotide exchange factor 7 |
| Target Name        | Rho guanine nucleotide exchange factor 7                             |
| Other Names        | Beta-Pix,PAK-interacting exchange factor beta,Pak3bp,PIXb            |
| Accession No.      | O55043   |

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.23ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Arhgef7

Note: This product is for in vitro research use only and is not intended for use in humans or animals.