

Mouse Cpm/Carboxypeptidase M ELISA Kit

Catalog No: #EK19651



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cpm/Carboxypeptidase M ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carboxypeptidase M
Target Name	Carboxypeptidase M
Other Names	CPM
Accession No.	Q80V42

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.098ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cpm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.