

Mouse Actn4/Alpha-actinin-4 ELISA Kit

Catalog No: #EK19657



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Actn4/Alpha-actinin-4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-actinin-4
Target Name	Alpha-actinin-4
Other Names	Non-muscle alpha-actinin 4
Accession No.	P57780

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.803ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Actn4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.