

# Rat Chi3l1/Chitinase-3-like protein 1 ELISA Kit

Catalog No: #EK19660



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Chi3l1/Chitinase-3-like protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chitinase-3-like protein 1
Target Name	Chitinase-3-like protein 1
Other Names	Cartilage glycoprotein 39,CGP-39
Accession No.	Q9WTV1

## Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.804ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Chi3l1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.