

Mouse Adam12/Disintegrin and metalloproteinase domain-containing protein 12 ELISA Kit

Catalog No: #EK19661



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Adam12/Disintegrin and metalloproteinase domain-containing protein 12 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 12
Target Name	Disintegrin and metalloproteinase domain-containing protein 12
Other Names	ADAM 12, Meltrin-alpha, Mltna
Accession No.	Q61824

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.47ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adam12

Note: This product is for in vitro research use only and is not intended for use in humans or animals.