

# Mouse Arid4b/AT-rich interactive domain-containing protein 4B ELISA Kit

Catalog No: #EK19662

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Arid4b/AT-rich interactive domain-containing protein 4B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse AT-rich interactive domain-containing protein 4B
Target Name	AT-rich interactive domain-containing protein 4B
Other Names	ARID domain-containing protein 4B, 180 kDa Sin3-associated polypeptide, Sin3-associated polypeptide p180, Histone deacetylase complex subunit SAP180, Sap180
Accession No.	A2CG63

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.186ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Arid4b

Note: This product is for in vitro research use only and is not intended for use in humans or animals.