

# Mouse Col10a1/Collagen alpha-1(X) chain ELISA Kit

Catalog No: #EK19664



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Col10a1/Collagen alpha-1(X) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(X) chain
Target Name	Collagen alpha-1(X) chain
Accession No.	Q05306

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Col10a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.