

Mouse Col7a1/Collagen alpha-1(VII) chain ELISA Kit

Catalog No: #EK19665



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Col7a1/Collagen alpha-1(VII) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(VII) chain
Target Name	Collagen alpha-1(VII) chain
Other Names	Long-chain collagen,LC collagen
Accession No.	Q63870

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.19ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Col7a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.