

Mouse Bmp15/Bone morphogenetic protein 15 ELISA Kit

Catalog No: #EK19666



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Bmp15/Bone morphogenetic protein 15 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein 15
Target Name	Bone morphogenetic protein 15
Other Names	BMP-15,Growth/differentiation factor 9B,GDF-9B,Gdf9b
Accession No.	Q9Z0L4

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.192ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bmp15

Note: This product is for in vitro research use only and is not intended for use in humans or animals.