

Human CD248/Endosialin ELISA Kit

Catalog No: #EK19667



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CD248/Endosialin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Endosialin
Target Name	Endosialin
Other Names	Tumor endothelial marker 1,CD248,CD164L1,TEM1
Accession No.	Q9HCU0

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.26pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CD248

Note: This product is for in vitro research use only and is not intended for use in humans or animals.