

Mouse Cdh2/Cadherin-2 ELISA Kit

Catalog No: #EK19668



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cdh2/Cadherin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cadherin-2
Target Name	Cadherin-2
Other Names	Neural cadherin, N-cadherin, CD325
Accession No.	P15116

Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.279ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdh2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.