

# Dog CELA1/Chymotrypsin-like elastase family member 1 ELISA Kit

Catalog No: #EK19672



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog CELA1/Chymotrypsin-like elastase family member 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Chymotrypsin-like elastase family member 1
Target Name	Chymotrypsin-like elastase family member 1
Other Names	Elastase-1,ELA1,ELS1
Accession No.	Q867B0

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/ml Sensitivity: 21pg/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CELA1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.