

Rat Cela1/Chymotrypsin-like elastase family member 1 ELISA Kit

Catalog No: #EK19674



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cela1/Chymotrypsin-like elastase family member 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chymotrypsin-like elastase family member 1
Target Name	Chymotrypsin-like elastase family member 1
Other Names	Elastase-1,Ela1
Accession No.	P00773

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:35pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cela1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.