

# Pig CALR/Calreticulin ELISA Kit

Catalog No: #EK19675



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig CALR/Calreticulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Calreticulin
Target Name	Calreticulin
Other Names	CRP55, Calregulin, Endoplasmic reticulum resident protein 60, ERp60, HACBP
Accession No.	P28491

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CALR

Note: This product is for in vitro research use only and is not intended for use in humans or animals.