

# Rat Adra2a/Alpha-2A adrenergic receptor ELISA Kit

Catalog No: #EK19679



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Adra2a/Alpha-2A adrenergic receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-2A adrenergic receptor
Target Name	Alpha-2A adrenergic receptor
Other Names	Alpha-2A adrenoreceptor,Alpha-2A adrenoceptor,Alpha-2D adrenergic receptor,CA2-47
Accession No.	P22909

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:8.02pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Adra2a

Note: This product is for in vitro research use only and is not intended for use in humans or animals.