

Mouse Ang2/Angiogenin-2 ELISA Kit

Catalog No: #EK19680



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ang2/Angiogenin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiogenin-2
Target Name	Angiogenin-2
Other Names	Angiogenin-related protein,Angrp
Accession No.	Q64438

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.17ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ang2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.