

Human C6/Complement component C6 ELISA Kit

Catalog No: #EK19682



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human C6/Complement component C6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Complement component C6
Target Name	Complement component C6
Accession No.	P13671

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C6

Note: This product is for in vitro research use only and is not intended for use in humans or animals.