

Mouse Bcl11a/B-cell lymphoma/leukemia 11A ELISA Kit

Catalog No: #EK19684

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Bcl11a/B-cell lymphoma/leukemia 11A ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse B-cell lymphoma/leukemia 11A |
| Target Name | B-cell lymphoma/leukemia 11A |
| Other Names | BCL-11A,B-cell CLL/lymphoma 11A,COUP-TF-interacting protein 1,Ecotropic viral integration site 9 protein,EVI-9,Ctip1,Evi9,Kiaa1809 |
| Accession No. | Q9QYE3 |

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.312ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bcl11a

Note: This product is for in vitro research use only and is not intended for use in humans or animals.