

# Mouse Bcl11b/B-cell lymphoma/leukemia 11B ELISA Kit

Catalog No: #EK19685

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bcl11b/B-cell lymphoma/leukemia 11B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse B-cell lymphoma/leukemia 11B
Target Name	B-cell lymphoma/leukemia 11B
Other Names	BCL-11B,B-cell CLL/lymphoma 11B,COUP-TF-interacting protein 2,Radiation-induced tumor suppressor gene 1 protein,mRit1,Ctip2,Rit1
Accession No.	Q99PV8

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.188ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bcl11b

Note: This product is for in vitro research use only and is not intended for use in humans or animals.