

Mouse Ca3/Carbonic anhydrase 3 ELISA Kit

Catalog No: #EK19686



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ca3/Carbonic anhydrase 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carbonic anhydrase 3
Target Name	Carbonic anhydrase 3
Other Names	Carbonate dehydratase III,Carbonic anhydrase III,CA-III,Car3
Accession No.	P16015

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.061ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ca3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.