

Rat Abca7/ATP-binding cassette sub-family A member 7 ELISA Kit

Catalog No: #EK19689



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Abca7/ATP-binding cassette sub-family A member 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat ATP-binding cassette sub-family A member 7
Target Name	ATP-binding cassette sub-family A member 7
Accession No.	Q7TNJ2

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Abca7

Note: This product is for in vitro research use only and is not intended for use in humans or animals.