

# Mouse Atp2a1/Sarcoplasmic/endoplasmic reticulum calcium ATPase 1 ELISA Kit

Catalog No: #EK19692



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Atp2a1/Sarcoplasmic/endoplasmic reticulum calcium ATPase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Target Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Other Names	SERCA1, Calcium pump 1, Calcium-transporting ATPase sarcoplasmic reticulum type, fast twitch skeletal muscle isoform, Endoplasmic reticulum class 1/2 Ca(2+) ATPase
Accession No.	Q8R429

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Atp2a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.