

# Human CASQ1/Calsequestrin-1 ELISA Kit

Catalog No: #EK19694



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CASQ1/Calsequestrin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calsequestrin-1
Target Name	Calsequestrin-1
Other Names	Calmitine, Calsequestrin, skeletal muscle isoform, CASQ
Accession No.	P31415

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.106ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CASQ1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.