

Mouse Casq1/Calsequestrin-1 ELISA Kit

Catalog No: #EK19695



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Casq1/Calsequestrin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calsequestrin-1
Target Name	Calsequestrin-1
Other Names	Calmitine, Calsequestrin, skeletal muscle isoform
Accession No.	O09165

Application Details

Assay Type: Sandwich Detect Range: 0.078-5ng/mL Sensitivity: 0.03ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Casq1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.