

# Rat Aldh4a1/Delta-1-pyrroline-5-carboxylate dehydrogenase, mitochondrial ELISA Kit

Catalog No: #EK19697

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Aldh4a1/Delta-1-pyrroline-5-carboxylate dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Delta-1-pyrroline-5-carboxylate dehydrogenase, mitochondrial
Target Name	Delta-1-pyrroline-5-carboxylate dehydrogenase, mitochondrial
Other Names	P5C dehydrogenase, Aldehyde dehydrogenase family 4 member A1, L-glutamate gamma-semialdehyde dehydrogenase
Accession No.	P0C2X9

## Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.37ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Aldh4a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.