

Rat Chad/Chondroadherin ELISA Kit

Catalog No: #EK19699



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Chad/Chondroadherin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chondroadherin
Target Name	Chondroadherin
Other Names	Cartilage leucine-rich protein
Accession No.	O70210

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.06ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Chad

Note: This product is for in vitro research use only and is not intended for use in humans or animals.