

Bovine COLEC11/Collectin-11 ELISA Kit

Catalog No: #EK19700



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Bovine COLEC11/Collectin-11 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Bovine |
| Specificity | Natural and recombinant Bovine Collectin-11 |
| Target Name | Collectin-11 |
| Other Names | Collectin kidney protein 1,CL-K1 |
| Accession No. | Q17QH6 |

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.065ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:COLEC11

Note: This product is for in vitro research use only and is not intended for use in humans or animals.