

Mouse Cdk5r1/Cyclin-dependent kinase 5 activator 1 ELISA Kit

Catalog No: #EK19706



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cdk5r1/Cyclin-dependent kinase 5 activator 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase 5 activator 1
Target Name	Cyclin-dependent kinase 5 activator 1
Other Names	CDK5 activator 1, Cyclin-dependent kinase 5 regulatory subunit 1, TPKII regulatory subunit, Cdk5r, Nck5a
Accession No.	P61809

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdk5r1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.