

Mouse Cflar/CASP8 and FADD-like apoptosis regulator ELISA Kit

Catalog No: #EK19707

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cflar/CASP8 and FADD-like apoptosis regulator ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CASP8 and FADD-like apoptosis regulator
Target Name	CASP8 and FADD-like apoptosis regulator
Other Names	Caspase homolog,CASH,Caspase-eight-related protein,Casper,Caspase-like apoptosis regulatory protein,CLARP,Cellular FLICE-like inhibitory protein,c-FLIP,FADD-like antiapoptotic molecule 1,FLAME-1,Inhibitor of FLICE,I-FLICE,MACH-related inducer of toxicity,MRIT,Usurpin,Cash
Accession No.	O35732

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.054ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cflar

Note: This product is for in vitro research use only and is not intended for use in humans or animals.