

Human ADTRP/Androgen-dependent TFPI-regulating protein ELISA Kit

Catalog No: #EK19708



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ADTRP/Androgen-dependent TFPI-regulating protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Androgen-dependent TFPI-regulating protein
Target Name	Androgen-dependent TFPI-regulating protein
Other Names	C6orf105
Accession No.	NP_116133.1/ncbi

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/ml Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ADTRP

Note: This product is for in vitro research use only and is not intended for use in humans or animals.