

Mouse Adcyap1r1/Pituitary adenylate cyclase-activating polypeptide type I receptor ELISA Kit

Catalog No: #EK19709

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Adcyap1r1/Pituitary adenylate cyclase-activating polypeptide type I receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Pituitary adenylate cyclase-activating polypeptide type I receptor
Target Name	Pituitary adenylate cyclase-activating polypeptide type I receptor
Other Names	PACAP type I receptor
Accession No.	P70205

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.7pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adcyap1r1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.