

# Mouse Casr/Extracellular calcium-sensing receptor ELISA Kit

Catalog No: #EK19712



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Casr/Extracellular calcium-sensing receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Extracellular calcium-sensing receptor
Target Name	Extracellular calcium-sensing receptor
Other Names	CaSR, Parathyroid cell calcium-sensing receptor, PCaR1, Gprc2a
Accession No.	Q9QY96

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.09ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Casr

Note: This product is for in vitro research use only and is not intended for use in humans or animals.